

INSTANT SURFACE TEST

DRUGSWAB™

Test Features:

- *Affordable*
- *Easy-to-Use*
- *Results in Minutes*

Drugs Detected:

- *Fentanyl*
- *Cocaine*
- *Opiates / Heroin*
- *Methamphetamine / Ecstasy*
- *Marijuana*



Testing Options:

Substances—Add the substance to the buffer vial. Shake for 15-30 seconds. Apply 15-20 drops from the buffer vial to the device collection pad

Drug Residue—Remove the blue cap and lightly brush the collection pad on the surface to collect the drug residue. Apply 15-20 drops from the buffer vial to the device collection pad.



800-668-0031

info@tryrapid.com

www.Tryrapid.com